

## NEW FEATURES 2015

### Advanced Training Features

SmartMan Pro Training continues to set the standard for CPR training. It now has those extra features that trainers have asked for and they are easy to access quickly.

All features are available immediately. Click on the button and they appear. Click on the button again to hide the feature. The trainer determines when the situation requires that feature to be displayed or not.

#### PAUSE ON/OFF

A single click to Pause an Activity. This allows you to explain a particular point, highlight a point or to correct a procedure.

A dark blue rectangular button with the word "Pause" in white text.


Click UnPause and you continue the activity where you left off.

For example, this is useful if you want to discuss a change over procedure during CPR. Or if you are doing a code recertification, you can pause to clarify roles and procedures for particular people without having to start again.

[More Information](#)

#### FEEDBACK ON/OFF

Simply Click the Button and the Real Time Color Feedback is Hidden from view. When the activity has finished, the score will be displayed and all of the colored bars will be displayed

A dark blue rectangular button with the text "Feedback On" in white.

Self Test: This is a great feature when a person is consistently achieving high scores. Simply click to turn off the feedback and the person can self-test to see how well they can do without the visual feedback and whether they have established muscle memory for correct performance.

Objective Evaluation: It can also be used as an objective evaluation of how well a person performs. The results record that the activity was done without the benefit of feedback.

[More Information](#)

## METRONOME ON/OFF



A Feature which you asked for in order to help those people who are having difficulty getting the rate correct.

Click on the Metronome Button. Then click Play. Stop it at any time by clicking on Stop.

The Metronome defaults to a rate of 120 per minute, but you can change the rate by using the slider.

[More Information](#)

## PARAMETERS

This feature provides maximum flexibility to the trainer. It is for trainers who wish to have students perform variations on the 30:2 compressions ratio.

Parameters

Set the number of compressions per cycle and the number of cycles.

For example set 1 cycle of 200 compressions. Or set 4 cycles of 50 compressions. SmartMan will calculate the times and the quality and will still provide students with an overall score and with analysis of the performance.

[More Information](#)

