

SmartMan Manual – 2. Actions Menu 6. Activate SmartMan

V5.1 and later

The Actions Menu

The Actions Menu is where you perform actions which are required by the program and you can view previous results of activities and tests performed by the person who is logged in.

Click on any item to perform the required action.

Actions Menu

1. Register
2. Login Next User
3. View Previous Results
4. Check for Updates
5. Calibrate SmartMan
6. Activate SmartMan
7. SmartMan ONLINE

6. Activate SmartMan

To use SmartMan, both the physical manikin and the software need to be linked.

You must activate the SmartMan program. To do this you must have a key code which is issued by SmartMan. If you do not have one, please contact us at: info@ambusmartman.com.

Please include the date of purchase, the company contacts, an email address for the key contact person and a telephone number.